

DL6MB Large Capacity Refrigerated Centrifuge(Blood Bank Centrifuge)



Features

1. Widely used in the field of blood bank, pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. There are 6 kinds of rotors for your choice.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

DL6MB Technical Parameter:

Max. Speed	6000rpm
Max. RCF	6600×g
Max. Capacity	6×1000ml(6*500ml blood bag)
Timer	1~9h59min
RPM/RCF Convert	Yes
Noise (dB)	≤ 60
Temperature Range	-20~40℃
Acc/Dec	10 Kinds
Speed Accuracy	±20r/min
Temperature Accuracy	±1℃
Voltage(V/Hz)	AC 220V 50HZ/60HZ
Size (L x W x Hmm)	800×740×930mm
Net Weight(Kg)	310KG
Certificates	CE,ISO & Calibration report are available

Matched Rotors for DL6MB

Order No.	Rotor	Max Speed (rpm)	Max Volume (ml)	Max.RCF (×g)
6MB-1	Swing Rotor	4200	6×1000ml	5100
6MB-2	Angle Rotor	6000	4×300ml	5390
6MB-3	Angle Rotor	6000	6×300ml	5660
6MB-4	Angle Rotor	6000	6×500ml	6600

