

TXK4 Gel Card (ID Card) Centrifuge



Application:

TXK4 blood group serological test multipurpose centrifuge is widely used in the test of blood group serology, blood type, microcolumn agglutination, immunodetection and etc.

Features:

- 1. Brushless motor to make sure free maintenance, no powder pollution, micro computer control and digital display the time, speed and RCF to let operation easily.
- 2. Special centrifugation technical to different reagent according to the requirement.
- 3. Regular centrifugation program, customers can press the button to choose the required one.
- 4. Electrical lid lock, over speed and imbalance protection.

TXK4 Technical Parameter:

Max. Speed	4000rpm		
Max. RCF	$2120 \times g(24Cards)$, $1480 \times g(12Cards)$		
Max. Capacity	24 Cards/12Cards		
Time Range	0~99min59sec		
RPM/RCF Convert	Yes		
Noise (dB)	≤ 55		
Temperature	Normal		
Acc/Dec	10 Kinds		
Speed Accuracy	±20r/min		
Voltage(V/Hz)	AC 220V/110V 50HZ/60HZ		
Size (W x D x Hmm) 12Cards	485×320×250mm		
Net Weight(Kg) 12 Cards	17KG		
Size (W x D x Hmm) 24 Cards	500×360×250mm		
Net Weight(Kg) 24 Cards	19KG		
Certificates	CE,ISO & Calibration report are available		

Matched Rotors for TXK4

Order No.	Rotor	Max Speed (rpm)	Max Volume (ml)	Max. RCF (×g)
K4-1	Swing Rotor	4000	12 cards	1480
K4-2	Swing Rotor	3000	24 cards	1190



