

GL10 High Speed Refrigerated Centrifuge



Features

1. Widely used in the field of pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. Kinds of rotors for your choice.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

GL10 Technical Parameter:

Max. Speed	10000rpm
Max. RCF	11900 × g
Max. Capacity	6 × 1000ml
Time Range	0~9h59min
RPM/RCF Convert	Yes
Noise (dB)	≤ 55
Temperature	-20~40℃
Acc/Dcc	10 Kinds
Speed Accuracy	± 20r/min
Temperature Accuracy	± 1℃
Voltage(V/Hz)	AC 380V/220V 50HZ/60HZ
Size (W x D x Hmm)	940 × 890 × 1000mm
Net Weight(Kg)	380KG
Certificates	CE,ISO & Calibration report are available

Matched Rotors for GL10

Order No.	Rotor No.	Max Speed (rpm)	Max Volume (ml)	Max.RCF (× g)
N030324	Angle Rotor	7000	6×1000ml	11900
N030123	Angle Rotor	8000	6×500ml	11740
N030122	Angle Rotor	10000	6×300ml	15730
N030298	Swing Rotor	4200	6×1000ml	5100